



PRODUCT BROCHURE · BENCHMARK SCIENTIFIC

A1700 (Agarose LE) Brochure

2 models inside · Labware & Glassware · Agarose



QUICK VIEW — SERIES AT A GLANCE

2 models in this document from Benchmark Scientific. Full specifications inside.

MODELS IN THIS DOCUMENT

A1700 · A1700-GROUP

SHIP WEIGHT

1 lb

SOURCING

Authorized U.S. distributor

Distributed by Legacy Lab Supplies

legacylabsupplies.com · (855) 746-7557 / 855-SHOP-LLS · quotes in 24 hours



AGAROSE LE

High Performance, General Purpose Agarose



Green Choice
Organic Solvent Free

Benchmark Agarose LE is a highly purified agarose, suitable for a variety of molecular biology applications. It is refined using an advanced process that *excludes the use of organic solvents*, yielding a cleaner end-product with a significantly reduced environmental impact. Agarose LE can be used for analyses of proteins and nucleic acids of various sizes (150 bp to 25 kbp).

It's low EEO ($\leq 0.13\text{-m}_r$) promotes increased electrophoretic mobility, yielding improved resolution and shorter run times. This also allows macromolecules and larger particles (subcellular fragments, viruses, etc.) to migrate more freely through the gel matrix. The consistently low EEO also provides a reduction in band distortion (caused by counterflow) that can result from the presence of excessive sulfate-rich negative ions.

Agarose LE is widely used for nucleic acid electrophoresis (analytical and preparative), protein electrophoresis (including radial immunodiffusion) and various blotting protocols. It is easily soluble, free of nucleases, and easy to use. It is highly transparent (forms a clear, colorless solution at 1g:100ml H_2O), and exhibits exceptionally low absorption of chemical staining agents. Pore size can be adjusted by simple modifications to the concentration ratio.

Formulated for high gel strength and integrity, Agarose LE exhibits exceptional thermal stability and mechanical resistance, ensuring safe, easy handling, regardless of whether a denaturing agent has been added.

- Multi-purpose, high purity
- Enhanced resolution and clarity
- Reduced background
- Low EEO/increased electrophoretic mobility
- RNase, DNase, Protease-free
- High hysteresis, improved thermal stability

Technical Data

Moisture:	$\leq 7\%$
EEO (electroendosmosis):	$\leq 0.13\text{-m}_r$
Sulfate:	$\leq 0.20\%$
Gelling Temp. (at 1.5%):	$36^\circ \pm 1.5^\circ \text{C}$
Melting Temp. (at 1.5%):	$88^\circ \pm 1.5^\circ \text{C}$
Gel Strength (at 1.0%):	1200g/cm^2
Gel Strength (at 1.5%):	2500g/cm^2
DNase, RNase:	None detected
Endonuclease:	None detected
Protease:	None detected

Ordering Information

A1700	Agarose LE, organic solvent free, 25g
A1701	Agarose LE, organic solvent free, 100g
A1705	Agarose LE, organic solvent free, 500g

Benchmark
Scientific

PO Box 709 Edison, NJ 08817
Phone: 908-769-5555 Fax: 732-313-7007
Web: www.BenchmarkScientific.com
Email: Info@BenchmarkScientific.com