



PRODUCT BROCHURE · BENCHMARK SCIENTIFIC

R3005.3010 (TubeRoller) Brochure

2 models inside · Mixing & Agitation · Rollers & Rotators



QUICK VIEW — SERIES AT A GLANCE

2 models in this document from Benchmark Scientific. Full specifications inside.

MODELS IN THIS DOCUMENT

R3005 · R3005-E

SHIP WEIGHT

12 lb

VOLTAGE

115 V / 60 Hz or 230 V / 50 Hz

SOURCING

Authorized U.S. distributor

Distributed by Legacy Lab Supplies

legacylabsupplies.com · (855) 746-7557 / 855-SHOP-LLS · quotes in 24 hours




TUBE ROLLER

TILTING ROLLER FOR TUBES AND BOTTLES



R3010

- + Gentle rolling motion
- + Smooth $\pm 5^\circ$ tilt action
- + For tubes/bottles from 8mm to 200mm diameter



R3005

Benchmark TubeRollers are designed for the gentle, thorough mixing of a variety of containers, including Vacutainers™, bottles, test tubes, centrifuge tubes and roller bottles. Both models (5 roller and 10 roller) operate at a constant 38rpm, chosen as the optimum speed for the broadest range of mixing requirements.


As the roller mechanism rotates through each complete revolution, the right side of each roller is gently elevated about 0.5 in (13mm) and returns to the horizontal starting position. This additional, gentle up/down motion provides a "wave" effect and guarantees that the entire sample is subjected to complete, thorough mixing, end to end. This built in tilt feature provides a level of homogeneous mixing uniformity that can not be matched by any stationary roller apparatus.

From blood collection tubes to roller bottles, Benchmark TubeRollers are the ideal choice for research, clinical laboratories and doctors' offices.


Ordering Information

R3005* TubeRoller with 5 rollers, 115V
R3010* TubeRoller with 10 rollers, 115V

* 115V with US plug. For 230V (EU plug), please add (-E)



R3005 shown with Vacutainers™



R3005 shown with 15ml conical tubes

Typical Sample Load Capacity:

Vessel	R3005	R3010
12 x 75mm tubes	12	26
16 x 100mm tubes	8	18
16 x 125mm tubes	8	18
15ml screw cap tubes	8	18
50ml screw cap tubes	8	18
850cm ² roller bottles	1	2

Technical Data

Speed: Fixed, 38 rpm ($\pm 10\%$)
 Motion: Simultaneous roll & tilt ($\pm 5^\circ$)
 Operating Temp: $+4^\circ$ to 40°C (external ambient temp.)

	R3005	R3010
Number of Rollers:	5	10
Dimensions (cm):	45 x 18 x 9	45 x 36 x 9
Dimensions (in):	18 x 7 x 3.5	18 x 14 x 3.5
Weight:	5kg / 11 lbs	8kg / 17.6 lbs

Individual Roller Dims.: 3 x 32 cm (1.2 x 12.5 in)
 Vessels Accepted: 8mm to 200mm diameter
 Roller Spacing: 5mm; rollers are on 35mm centers
 Electrical: 120V, 60Hz or 230V, 50-60Hz
 Warranty: 2 years

PO Box 709
Edison, NJ 08818
Phone: 908-769-5555
Fax: 908-222-1864
Web: www.BenchmarkScientific.com
Email: Info@BenchmarkScientific.com

Benchmark 
 Scientific